



PHILIPS
RESPIRONICS

DreamStation

Auto BiPAP

Bi-level therapy for OSA patients: a comfortable alternative to CPAP

DreamStation Auto BiPAP features

- Advanced event detection*
- Therapy modes: CPAP, fixed bi-level, auto bi-level
- Adjustable PS min and PS max
- **NEW** - Smart Ramp (0-45 min)
- Flex Comfort: C-Flex, C-Flex+, Bi-Flex
- **NEW** - Performance Check
- Encore patient management system: compliance and data reports, waveforms reports
- **NEW** - Multiple remote monitoring solutions including integrated Bluetooth® technology and optional cellular or wifi modems
- On-board memory
- Oximetry module (optional)

Humidification

- Active humidification controlled by the room environment
- Pre-heat humidifier capable
- Heated tube (optional)

DreamStation Auto BiPAP benefits

- **Daily Progress Feedback** displays the length of sleep and progress, designed to motivate the patient in the early days of therapy.
- **Performance Check** simplifies in-home device evaluation for easy remote troubleshooting.
- **Connectivity:** DreamStation offers more remote monitoring options than any leading CPAP brand.

*Philips Respiromics devices with advanced event detection detect, and log these advanced events: Open airway apnea (CA), Obstructed airway apnea (OA), Hypopnea (H), respiratory events-related arousal (RERA), flow limitation (FL), vibratory snore (VS), large leak (LL), periodic breathing (PB).

DreamStation Auto BiPAP specifications

Part number	Description
INX700S15	DreamStation Auto BiPAP with SD card
INX700H15	DreamStation Auto BiPAP with humidifier and SD card

Specifications

Therapy modes	Fixed CPAP, fixed BiPAP, Auto BiPAP
Therapy features	EZ-Start
Flex pressure relief	0 to 3
Pressure range	4-25 cm H ₂ O
Smart Ramp time	0 to 45 min (5-min increments)
Starting ramp pressure	Patient adjustable
Humidification	Heated humidification: fixed, adaptive, heated tube modes Heated tube (optional)
Device controls	Color LCD/control wheel/push button
LCD backlight control	Ambient light sensor
Filters	Reusable pollen and disposable ultra-fine
Data storage capacity (minimum)	Display: compliance 1 day, 7-day, 30-day averages SD card > 1 year On-board: up to 6 months
Altitude compensation	Automatic
Warranty	2 years
Dimensions	15,7 (L) x 19,3 (W) x 8,4 cm (H) (without humidifier) 29,7 (L) x 19,3 (W) x 8,4 cm (H) (with humidifier)
Weight	1,33 kg (without humidifier), 1,98 kg (with humidifier) but including power supply
Power supply	80W/12V/6,67A
Electrical requirements	100-240 VAC, 50/60 Hz, 2.0-1.0 A

Sleep therapy accessories

Humidification	Part number
DreamStation heated humidifier	EUXH
Water chamber	1122520

Tubing

Tubing	Part number
Heated tube (15 mm)	HT15
Standard tube (15 mm)	PR15

Filters

Filters	Part number
Pollen (foam) filter, reusable (1/pack)	1122446
Ultra-fine filter, disposable (1/pack)	1122447
Ultra-fine filter, disposable (2/pack)	1122518
Ultra-fine filter, disposable (6/pack)	1122519

Accessories

Accessories	Part number
DreamStation carrying case	1121162
DreamStation link module/SpO ₂	1120293
DreamStation SpO ₂ cable	1121694
Reuseable finger sensor (clip)	936
Oximetry finger sensor, Adult, Flex	953
Flexiwrap, sens tape, Adult (25/pack)	954A
Reuseable finger sensor (soft), medium	1062050
Reuseable finger sensor (soft), small	1062051
DreamStation, 80W power supply	1118499

Connectivity

Connectivity	Part number
DreamStation cellular/SpO ₂ modem	100621C
DreamStation Wifi/SpO ₂ modem	100721W





PHILIPS

RESPIRONICS

DreamStation

Auto CPAP

Designed to **empower patients** and facilitate APAP therapy adoption

DreamStation Auto CPAP features

- Advanced event detection*
- Clinically validated auto algorithm (responds to obstructive events: AO, H, RE, FL, VS.)
- Proactive search algorithm
- CPAP-Check
- Auto-Trial
- Opti-Start
- **NEW** - EZ-Start
- **NEW** - Smart Ramp (0-45 min)
- A-Flex
- **NEW** - Performance Check
- Encore patient management system: compliance and data reports, waveforms reports
- **NEW** - Multiple remote monitoring solutions including cellular, Bluetooth® and Wifi modems
- On-board memory
- Oximetry module (optional)

Humidification

- Active humidification controlled by the room environment
- Pre-heat humidifier capable
- Heated tube (optional)

DreamStation Auto CPAP benefits

- **EZ-Start's** automatic, personalized adjustments to CPAP pressure is designed to help patients gradually acclimate to their prescribed level of therapy.
- **Daily Progress Feedback** displays the length of sleep and progress, designed to motivate the patient in the early days of therapy.
- **Performance Check** simplifies in-home device evaluation for easy remote troubleshooting.
- **Connectivity:** DreamStation offers more remote monitoring options than any leading CPAP brand.

*Philips Respironics devices with advanced event detection detect, and log these advanced events: Open airway apnea (CA), Obstructed airway apnea (OA), Hypopnea (H), respiratory events-related arousal (RERA), flow limitation (FL), vibratory snore (VS), large leak (LL), periodic breathing (PB).

DreamStation Auto CPAP specifications

Part number	Description
INX500S15	DreamStation Auto CPAP with SD card
INX500H15	DreamStation Auto CPAP with humidifier and SD card

Specifications

Therapy modes	Auto, CPAP, CPAP-Check
Therapy features	Auto-Trial, EZ-Start, Opti-Start
Flex pressure relief	0 to 3
Pressure range	4-20 cm H ₂ O
Smart Ramp time	0 to 45 min (5-min increments)
Starting ramp pressure	Patient adjustable
Humidification	Heated humidification: fixed, adaptive, heated tube modes Heated tube (optional)
Device controls	Color LCD/control wheel/push button
LCD backlight control	Ambient light sensor
Filters	Reusable pollen and disposable ultra-fine
Data storage capacity (minimum)	Display: compliance 1 day, 7-day, 30-day averages SD card > 1 year On-board: up to 6 months
Altitude compensation	Automatic
Warranty	2 years
Dimensions	15,7 (L) x 19,3 (W) x 8,4 cm (H) (without humidifier) 29,7 (L) x 19,3 (W) x 8,4 cm (H) (with humidifier)
Weight	1,33 kg (without humidifier), 1,98kg (with humidifier) but including power supply
Power supply	80W/12V/6,67A
Electrical requirements	100-240 VAC, 50/60 Hz, 2.0-1.0 A

Sleep therapy accessories

Humidification	Part number
DreamStation heated humidifier	EUXH
Water chamber	1122520

Tubing

Tubing	Part number
Heated tube (15mm)	HT15
Standard tube (15mm)	PR15

Filters

Filters	Part number
Pollen (foam) filter, reusable (1/pack)	1122446
Ultra-fine filter, disposable (1/pack)	1122447
Ultra-fine filter, disposable (2/pack)	1122518
Ultra-fine filter, disposable (6/pack)	1122519

Accessories

Accessories	Part number
DreamStation carrying case	1121162
DreamStation link module	1120293
DreamStation SpO ₂ module	1121694
Reuseable finger sensor (clip)	936
Oximetry finger sensor, Adult, Flex	953
Flexiwrap, sens tape, Adult (25/pack)	954A
Reuseable finger sensor (soft), medium	1062050
Reuseable finger sensor (soft), small	1062051
DreamStation, 80W power supply	1118499

Connectivity

Connectivity	Part number
DreamStation cellular/SpO ₂ modem	100621C
DreamStation Wifi/SpO ₂ modem	100721W





PHILIPS
RESPIRONICS

DreamStation

BiPAP Pro

DreamStation helps support long-term patient use while creating **new efficiencies for your business.**

Connectivity options offer more monitoring choices than any leading CPAP brand with integrated Bluetooth® on every device, and offering optional Wi-Fi and cellular modems.*

Performance Check simplifies in-home device evaluation and can help reduce the number of normally functioning devices that are returned for service.

EZ-Start's automatic, personalized adjustments to CPAP pressure can help patients gradually acclimate to their prescribed level of therapy.

Daily Progress Feedback displays a simple trend of patients' nightly use to help track progress on the path to adherence.

DreamStation features

- Advanced event detection
- Mask type resistance control
- Mask fit check
- Oximetry capable
- Encore patient management system with optional Adherence Profiler
- DreamMapper patient self-management system compatibility

Humidification

- Fixed, adaptive, heated tube
- Pre-heat humidifier option
- Modular design with one-piece, dishwasher-safe water chamber

Protection programs

- Flex Promise program
- Bi-level Rescue program
- "One-Touch" warranty programs

BiPAP Pro features

- C-Flex, C-Flex+, Bi-Flex
- SmartRamp
- EZ-Start
- CPAP-Check
- Advanced event detection
- Waveform data
- Standard integrated Bluetooth

*Internal assessment of 2015 competitive CPAP data comparing to ResMed Aircurve10/Aircurve10 platform and Fisher & Paykel Icon series platform. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and the use of such marks by Philips Respironics is under license. Other trademarks and trade names are those of their respective owners.

DreamStation BiPAP Pro

Comfort features			Encore data capture and reporting	
Modes	Features	Flex Comfort	Advanced Event Detection	Flow Waveforms
Fixed CPAP Fixed BiPAP	EZ-Start	C-Flex, C-Flex+, Bi-Flex	Standard	Standard

Specifications

DreamStation BiPAP Pro	
Pressure range	4 to 25 cm H2O
Flex pressure relief	0 to 3
Ramp time	0 to 45 min (5-minute increments)
Starting ramp pressure	4 cm H2O to CPAP Pressure (CPAP); 4cm H2O to EPAP (BiPAP)
Humidification	Heated humidification: Fixed, Adaptive, Heated Tube modes; Heated Tube (optional)
Data storage capacity (minimum)	SD card: > 1 year On-board: <6 months
Altitude compensation	Automatic
Filters	Reusable pollen; disposable ultra-fine
Device controls	LCD, control dial/push button
Warranty	2 years (US)
Dimensions	15.7x19.3x8.4cm (w/o humidifier) 29.7x19.3x8.4cm (w/ humidifier)
Weight	1.33kg/2.94lbs (w/o humidifier) 1.98kg/4.37lbs (w/humidifier) Incl. power supply
Electrical requirements	100 – 240 VAC, 50/60 Hz, 2.0-1.0 A

Ordering information

DreamStation BiPAP Pro	Part number
DreamStation BiPAP Pro	INX600S15
DreamStation BiPAP Pro with humidifier	INX600H15
DreamStation BiPAP Pro with humidifier and heated tube	INX600T15
Power supply	Part number
DreamStation 80-watt power supply	1118499
Power cord	1005894

Ordering information

Humidification	Part number
DreamStation humidifier	DSXH
DreamStation water chamber	1122520
RP-DS Humidifier Dry Box Assembly	1120668
RP-DS Humidifier Flip Lid Seal	1120617
RP-DS Humidifier Dry Box Inlet Seal	1120613
Tubes	Part number
15MM Std Tube – DreamStation - RP	PR15
15MM Heated Tube	HT15
Filters	Part number
Pollen filter, reusable (1 per pack)	1122446
Ultra-fine filter, disposable (1 per pack)	1122447
Ultra-fine filter, disposable (2 per pack)	1122518
Ultra-fine filter, disposable (6 per pack)	1122519
Accessories	Part number
Link module	1120293
Nonin SpO2, assembly	1121694
Reuseable finger sensor (clip)	936
Oximetry finger sensor, adult, flex	953
Flexiwrap, sens tape, adult, 25/pk	954A
Reuseable finger sensor (soft), medium	1062050
Reuseable finger sensor (soft), small	1062051
Modems	Part number
Cellular modem, oximetry capable	100600C
Cellular modem, non-oximetry	100605C
Wi-Fi accessory, oximetry capable	100700W
Wi-Fi accessory, non-oximetry	100705W





PHILIPS

RESPIRONICS

DreamStation

CPAP Pro

Designed to **empower patients** and facilitate CPAP therapy adoption

DreamStation CPAP Pro features

- Advanced event detection*
- Therapy modes: CPAP, CPAP-Check
- **NEW** - EZ-Start
- **NEW** - Smart Ramp (0-45 min)
- C-Flex
- **NEW** - Performance Check
- Encore patient management system: compliance and data reports, waveforms reports
- **NEW** - Multiple remote monitoring solutions including cellular, Bluetooth® and Wifi modems
- On-board memory
- Oximetry module (optional)

Humidification

- Active humidification controlled by the room environment
- Pre-heat humidifier capable
- Heated tube optional

DreamStation CPAP Pro benefits

- **EZ-Start's** automatic, personalized adjustments to CPAP pressure is designed to help patients gradually acclimate to their prescribed level of therapy.
- **Daily Progress Feedback** displays the length of sleep and progress, designed to motivate the patient in the early days of therapy.
- **Performance Check** simplifies in-home device evaluation for easy remote troubleshooting.
- **Connectivity:** DreamStation offers more remote monitoring options than any leading CPAP brand.

*Philips Respironics devices with advanced event detection detect, and log these advanced events: Open airway apnea (CA), Obstructed airway apnea (OA), Hypopnea (H), respiratory events-related arousal (RERA), flow limitation (FL), vibratory snore (VS), large leak (LL), periodic breathing (PB).

DreamStation CPAP Pro specifications

Part number	Description
INX400S15	DreamStation CPAP Pro with SD card
INX400H15	DreamStation CPAP Pro with humidifier and SD card

Specifications

Therapy modes	CPAP, CPAP-Check
Therapy features	EZ-Start
Flex pressure relief	0 to 3
Pressure range	4-20 cm H ₂ O
Smart Ramp time	0 to 45 min (5-min increments)
Starting ramp pressure	Patient adjustable
Humidification	Heated humidification: fixed, adaptive, heated tube modes Heated tube (optional)
Device controls	Color LCD/control wheel/push button
LCD backlight control	Ambient light sensor
Filters	Reusable pollen and disposable ultra-fine
Data storage capacity (minimum)	Display: compliance 1 day, 7-day, 30-day averages SD card > 1 year On-board: up to 6 months
Altitude compensation	Automatic
Warranty	2 years
Dimensions	15,7 (L) x 19,3 (W) x 8,4 cm (H) (without humidifier) 29,7 (L) x 19,3 (W) x 8,4 cm (H) (with humidifier)
Weight	1,33 kg (without humidifier), 1,98kg (with humidifier) but including power supply
Power supply	80W/12V/6,67A
Electrical requirements	100-240 VAC, 50/60 Hz, 2.0-1.0 A

Sleep therapy accessories

Humidification	Part number
DreamStation heated humidifier	EUXH
Water chamber	1122520

Tubing

Tubing	Part number
Heated tube (15mm)	HT15
Standard tube (15mm)	PR15

Filters

Filters	Part number
Pollen (foam) filter, reusable (1/pack)	1122446
Ultra-fine filter, disposable (1/pack)	1122447
Ultra-fine filter, disposable (2/pack)	1122518
Ultra-fine filter, disposable (6/pack)	1122519

Accessories

Accessories	Part number
DreamStation carrying case	1121162
DreamStation link module	1120293
DreamStation SpO ₂ module	1121694
Reuseable finger sensor (clip)	936
Oximetry finger sensor, Adult, Flex	953
Flexiwrap, sens tape, Adult (25/pack)	954A
Reuseable finger sensor (soft), medium	1062050
Reuseable finger sensor (soft), small	1062051
DreamStation, 80W power supply	1118499

Connectivity

Connectivity	Part number
DreamStation cellular/SpO ₂ modem	100621C
DreamStation Wifi/SpO ₂ modem	100721W

