

Description for Kiran Maxin fabric Radiation protection Lead Apron

Kiran Maxin fabric Radiation protection Lead Apron – the perfect combination of comfort and safety.

Made with lightweight, soft Maxin fabric, this lead apron offers maximum radiation protection while reducing fatigue during long procedures.

The apron is waterproof and available in three attractive colors – Azure blue, Graphite grey and Ocean green.

With a Lead Equivalence ranging from 0.25 to 0.50mm, you can feel safe and comfortable wearing this convenient radiation shielding gear.

Features of Kiran Maxin fabric Radiation protection Lead Apron :

- Comfortable
- Light weight
- Feather soft
- Maxin fabric
- Easy usage
- Reduced fatigue
- Available in 3 colours - Azure blue, Graphite grey and Ocean green

Specifications of Kiran Maxin fabric Radiation protection Lead Apron :

Product Name	Kiran Maxin fabric Radiation protection Lead Apron
Brand	Kiran
Model	Radiation protection Lead Apron
Pattern	Plain
Sizes Available	Medium, Large, Free Size, Small
Lead Equivalence	0.25 mm and 0.50 mm
Material	Maxin fabric
Color	Azure blue, graphite grey and ocean green
Application	Hospital, Laboratory
Feature	Waterproof